


EVRL

Erchonia Violet and Red Laser

Pivoting Laser Head to Aid Device Application

40 Degrees of Laser Head Movement

Easy to Use Touchscreen Technology



Only 7" Tall

Sleek Design Less than 1" Thick



Introducing the Newest Line of Cold Lasers from Erchonia®

The EVRL Violet Laser System is the newest Erchonia handheld laser in "smart phone" size with a user-friendly graphical user interface (GUI) and a sleek new cord free design. Operating at a wavelength of 405nm and 635nm, the EVRL Laser provides the practitioner with an innovative tool for treatment versatility.

The EVRL Laser System Advantage:

- **Anti-bacterial & Anti-microbial**
- **Red & Violet Laser Design**
- **Touchscreen GUI**

SPECIFICATIONS

Configuration: (1) 7.5mW Line Generated Red Laser Diode, (1) <5mW Line Generated Violet Laser Diode with Patented Optics

Wavelengths: 635nm & 405nm

Modulation: Constant Wave (CW) & Pulsed (Hz) (Red Only)

40 User Defined Memory Channels

Display: Full Color TFT Touchscreen Technology for Ease of Use

Power Source: 3.7 VDC Lithium-ion Polymer Battery

Enclosure: Micro-Textured Anodized Aircraft Aluminum for Strength and Ease of Cleaning

Weight: 0.95lbs (.42 kg) (Device Only)

Accessories: Power Supply, Wireless Charging Base, Carrying Case and Laser Safety Glasses

Compliant to ISO 13485 Medical Device Quality, IEC 60601-1 Safety, 60601-1-2 EMC, IEC 60825-1 Laser Safety, CE Mark

Laser Class 2 / Device Class II (USA); 2a (EU)

US PAT 6,013,096; US PAT 6,746,473 For additional US and International patents and patent pending information go to www.erchonia.com.



ERCHONIA®
World Leaders in Low-Level Laser Technology

(888)242-0571 Toll Free

www.erchonia.com

